



LAP-iX
Suction and
Irrigation Electrode



ADVANCING SURGICAL POSSIBILITIES



LAP-ix
Suction Irrigation Electrode

Electrode

S/I Only



L-Hook Electrode Tip

The sharp contact surface makes easy electric cutting. Rounded non-contact surface prevents unintended damage during pushing or lifting the device during the operation



Spoon Electrode Tip

Rounded shaped edge prevents unintended damage during the operation



Suction And Irrigation Tip

The metal tip is applied for the suction and irrigation tip to enhance the intensity

| TYPE | PRODUCT NO. | ELECTRODE | QUANTITY |
|---------|-------------|-----------|----------|
| Pistol | SI-PU-01-01 | L-Hook | 1ea/Box |
| | SI-PU-02-01 | Spoon | 1ea/Box |
| | SI-PO | S/I Only | 2ea/Box |
| Trumpet | SI-TU-01-01 | L-Hook | 1ea/Box |
| | SI-TU-02-01 | Spoon | 1ea/Box |
| | SI-TO | S/I Only | 2ea/Box |
| Tip | SI-ED-01-01 | L-Hook | 1ea/Box |
| | SI-ED-02-01 | Spoon | 1ea/Box |

Advancing surgical possibilities

LAWMED Ltd.

Unit 1, Walton Lodge, Bridge Street
Walton On Thames, KT12 1BT
United Kingdom

E: enquiries@lawmed.co.uk
T: +44 (0)1932 260 838
www.lawmed.co.uk

Manufactured by
SEJONG
MEDICAL

Distributed by
LAWMED Ltd.