



MID-TUBE™

OROGASTRIC CALIBRATION TUBE



36 Fr



38 Fr



48 Fr



MID™
INNOVATION FOR SURGERY

MID-TUBE™ FEATURES:

MID-TUBE™ is a silicone tube with a diameter of 36 Fr, 38 Fr or 48 Fr.

MID-TUBE™ can also be connected to a suction system for easy removal of air or fluid in the patient's stomach either before or during the procedure (leak test).

6 graduations (from 35 to 50 mm).

EXAMPLES OF PROCEDURE IN OBESITY SURGERY:

Gastric by-pass
Sleeve gastrectomy
Biliopancreatic diversion

...



TECHNICAL DATA:

	MID136	MID130	MID131
Material used	Medical Silicone		
Tube length	800 mm		791 mm
Tube external diameter	12 mm (36 Fr)	12.8 mm (38 Fr)	16 mm (48 Fr)

PACKAGING & PRODUCT CODE:

One unit packaged in a sealed pouch per box.

REF **MID130**

Shelf box dimensions: 24.6 x 17.5 x 3.4 cm.

REF **MID136**

Shelf box dimensions: 35 x 17.5 x 3.4 cm.

REF **MID131**

WARNINGS:



The MID-TUBE™ must only be handled by practitioners with experience of digestive surgery.



The particular instructions for use are shown in the notice.



Single Use / Do not resterilize.



Contents are sterile if package is unopened and undamaged.



Sterile (Ethylene Oxide).



Store in a cool, dark, dry place.

CE 0459 - Classe : IIa

MID 
INNOVATION FOR SURGERY

YOUR

LAWMED

ADVANCING SURGICAL POSSIBILITIES TOGETHER

NEED MORE INFORMATION
enquiries@lawmed.co.uk
+44 (0)1932 260 838

MEDICAL INNOVATION
DEVELOPPEMENT

9 CHEMIN DU JUBIN
69570 DARDILLY
FRANCE

Tel.: +33 (0)4 78 17 48 04

Fax: +33 (0)4 72 82 91 23

mid@mid-med.com

.....



www.mid-med.com