

## Disposable LESS Port

Move freely, smoothly and comfortably

- Instruments can be inserted under direct vision due to transparent platform
- Specially designed working channels help to avoid interaction of different instruments
- The instruments can be inserted freely and smoothly.
- Rotatable multi-port platform can provide different view angle for better manipulations.



Quick smoke evacuation,  
clear surgery field

The smoke evacuation tube length  
and position can be adjustable,  
it helps to increase the smoke  
evacuation efficiency and clear the  
surgery field

Single person installation,  
easy and reliable

- Installation can be done by single person;
- Buckle design makes the connection more reliable and secured.

## Transumbilical laparoendoscopic single site surgery

Product	Dimension	Code	Description
	Φ60	191-JQST4C635	glove type-5+5+10+12 -φ60*150
	Φ80	191-JQST4C835	glove type-5+5+10+12 -φ80*150
	Φ80	191-JQGJ4C835	silicone type-5+5+10+12 -φ80*150

## Vaginal natural orifice transluminal endoscopic surgery

Product	Dimension	Code	Description
	Φ60	191-JYST4C635	glove type-5+5+10+12 -φ60*150
	Φ80	191-JYST4C835	glove type-5+5+10+12 -φ80*150
	Φ80	191-JYGJ4C835	silicone type-5+5+10+12 -φ80*150

## Transanal minimally invasive surgery

Product	Dimension	Code	Description
	Φ80	191-JGST4C825	glove type-5+5+10+12 -φ80*90
	Φ80	191-JGST4C815	glove type-5+5+10+12 -φ80*70
	Φ80	191-JGGJ4C825	silicone type-5+5+10+12 -φ80*90
	Φ80	191-JGGJ4C815	silicone type-5+5+10+12 -φ80*70

## Disposable Wound Protectors/Retractors



### Protection:

- Decrease Surgical Site Infection Rate
- Prevent from cancer cell self-seeding
- Keep incision with proper moisture and temperature
- Help to stop bleeding

### Exposure:

- Maximize the port with the smallest incision size
- Better exposure without additional injury caused by metal retractors
- Reduce pressure and fatigue of using the traditional retractors

### Convenience:

- Easy Installation and dismantling
- Suitable for different surgeries



Protectors



Cap



Simple LESS Port

## Variable Height

Product	Dimension	Code	Incision sizes(mm)
	60/70-60/150	191-B06070645	20-40
	80/90-80/150	191-B08090845	25-60
	120/130-120/250	191-B12131255	50-90
	150/160-150/250	191-B15161555	70-120
	180/190-180/250	191-B18191855	90-140
	220/230-220/250	191-B22232255	110-170
	270/280-270/250	191-B27282755	140-220
	320/330-320/250	191-B32333255	170-250

## Variable Height With Cap (Code unchanged, sales code differentiated)

Product	Dimension	Code	Incision sizes(mm)
	60/70-60/150	191-B06070645	20-40
	80/90-80/150	191-B08090845	25-60
	120/130-120/250	191-B12131255	50-90
	150/160-150/250	191-B15161555	70-120

## Fixed Height

Product	Dimension	Code	Incision sizes(mm)
	40/30-15/25	191-D04030125	12-15
	50/40-25/25	191-D05040225	20-25
	70/70-35/40	191-D07070335	30-35
	60/50-25/5	191-D06050215	20-25
	100/50-40/5	191-D10050415	30-40