

Polymer Ligation Clip

Secure Ligation.
Enhanced Control.



Optimised Clip Stability
Enhanced clip security within the applicator jaws for improved control

Locking Mechanism
Secure closure with tactile and audible feedback to confirm proper engagement

Bow-Shaped Profile
Designed to allow safe removal using a compatible clip remover

Anti-Slip Ridges
Reliable 360° grip designed to minimise dislodgement and slippage

Advanced Hinge Design
Larger distal opening compared with conventional clips, enabling effective ligation through a trocar

NEED MORE INFORMATION

Please contact us at
enquiries@lawmed.co.uk or **+44 (0)1932 260 838**.

Lawmed Ltd
Building 5 The Heights
Wellington Way
KT13 0NY
UK


www.lawmed.co.uk







mindray

Manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

LM-S0015
© 2026 Lawmed

Polymer Ligation Clip



Product	Product Code	Dimension	UOM
	191-01636	ML, 6 Clips	10
	221-55636	5mm x 330	1
	191-01646	L, 6 Clips	10
	221-05646	10mm x 330	1
	191-01656	XL, 6 Clips	10
	221-05656	10mm x 330	1

NEED MORE INFORMATION

Please contact us at
enquiries@lawmed.co.uk or **+44 (0)1932 260 838**.

Lawmed Ltd
Building 5 The Heights
Wellington Way
KT13 0NY
UK

www.lawmed.co.uk

mindray

Manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

LM-S0015
© 2026 Lawmed