

Laparoscopic Retrieval Pouch



Features a reinforced seam and a flexible pouch, enabling the safe and efficient removal of tissue from the abdominal cavity during laparoscopic procedures.

It incorporates an automatic deployment mechanism and is available in a wide range of sizes to accommodate various tissue types.

All Prima bags are compatible with a 10 mm trocar, except for the XX-large size, which is introduced via a 12 mm trocar, eliminating the need to open an additional 15 mm trocar.

Product Codes	Product Description	Dimensions (L x W)	Volume	Trocar Size	Units Per Box
B21-1002	U- TYPE, Small	100mm x 100mm	200 ml	10 mm	5
B21-1003	U- TYPE, Medium	125mm x 125mm	410 ml	10 mm	5
B21-1004	U- TYPE, Large	160mm x 160mm	920 ml	10 mm	5
B21-1005	U- TYPE, X-Large	200mm x 200mm	1540 ml	10 mm	5
B11-1208	U- TYPE, XX-Large	300mm x 250mm	4800 ml	12 mm	5

NEED MORE INFORMATION

Please contact us at
enquiries@lawmed.co.uk or **+44 (0)1932 260 838**.

Lawmed Ltd
Building 5 The Heights
Wellington Way
KT13 0NY
UK

www.lawmed.co.uk