


# LAP SINGLE VISION

Laparoscopic Single Port



LS-51-1313      LS-51-1316      LS-51-1318

TYPE	PRODUCT NO.	DIMENSION	WOUND DIAMETER	QUANTITY
Normal	LS-51-1313	10mm x 3	60mm x 60mm	10ea/Box
	LS-51-1316	10mm x 3, 12mm x 1	60mm x 60mm	10ea/Box
	LS-51-1318	10mm x 3, 15mm x 1	60mm x 60mm	10ea/Box
	LS-51-3313	10mm x 3	60mm x 90mm	10ea/Box
	LS-51-3316	10mm x 3, 12mm x 1	60mm x 90mm	10ea/Box
	LS-51-3318	10mm x 3, 15mm x 1	60mm x 90mm	10ea/Box
Tube	LS-52-1313	10mm x 3	60mm x 60mm	10ea/Box
	LS-52-1316	10mm x 3, 12mm x 1	60mm x 60mm	10ea/Box
	LS-52-1318	10mm x 3, 15mm x 1	60mm x 60mm	10ea/Box
	LS-52-3313	10mm x 3	60mm x 90mm	10ea/Box
	LS-52-3316	10mm x 3, 12mm x 1	60mm x 90mm	10ea/Box
	LS-52-3318	10mm x 3, 15mm x 1	60mm x 90mm	10ea/Box



PRODUCT NO.	DIMENSION
SJ-SF-03	80mm x 40mm x 23mm, valve 270mm

**Reduce Odor And Mutagenic Chemicals**

Smoke, created in the burning of tissue contains mutagenic chemicals. Sejong Smoke Filter contains PRE/ACF/DNA Filters which reduces the chemicals and odor in the smoke

**Transparency**

Transparent material makes easy to access the abdominal cavity

**Quickly Removes The Smoke**

Isolation of gas inlet and smoke exhaust for smooth emission of the smoke in the abdominal

**Minimized Instrument Collision**

Minimizes the collision of the instrument during the single incision surgery, using silicon for the body



# LAPSINGLE PLUS

Laparoscopic Single Port




Lapsingle-41-1313 Lapsingle-41-1316 Lapsingle-41-1318

TYPE	PRODUCT NO.	DIMENSION	WOUND DIAMETER	QUANTITY
Normal	Lapsingle-41-1313	10mm x 3	60mm x 60mm	10ea/Box
	Lapsingle-41-1316	10mm x 3, 12mm x 1	60mm x 60mm	10ea/Box
	Lapsingle-41-1318	10mm x 3, 15mm x 1	60mm x 60mm	10ea/Box
Tube	Lapsingle-41-3316	10mm x 3, 12mm x 1	60mm x 90mm	10ea/Box
	Lapsingle-42-1313	10mm x 3	60mm x 60mm	10ea/Box
	Lapsingle-42-1316	10mm x 3, 12mm x 1	60mm x 60mm	10ea/Box
	Lapsingle-42-1318	10mm x 3, 15mm x 1	60mm x 60mm	10ea/Box
	Lapsingle-42-3316	10mm x 3, 12mm x 1	60mm x 90mm	10ea/Box



# LAPSINGLE

Laparoscopic Single Port



Lapsingle-31-1334 Lapsingle-31-1345 Lapsingle-31-1304 Lapsingle-31-1333

TYPE	PRODUCT NO.	DIMENSION	WOUND DIAMETER	QUANTITY
Normal	Lapsingle-31-1334	12mm x 1, 5mm x 3	60mm x 60mm	10ea/Box
	Lapsingle-31-1345	12mm x 2, 5mm x 2	60mm x 60mm	10ea/Box
	Lapsingle-31-1304	5mm x 4	60mm x 60mm	10ea/Box
	Lapsingle-31-1333	12mm x 1, 5mm x 2	60mm x 60mm	10ea/Box



# SEARCH PORT

Laparoscopic Single Port



TYPE	PRODUCT NO.	WOUND DIAMETER	QUANTITY
String	WR-21-6060	60mm x 60mm	10ea/Box
	WR-21-9090	90mm x 90mm	10ea/Box
	WR-21-9090V	60mm x 90mm	10ea/Box
Non-String	WR-22-6060	60mm x 60mm	10ea/Box
	WR-22-9090	90mm x 90mm	10ea/Box
	WR-22-9090V	60mm x 90mm	10ea/Box

Transparency  
Transparent material makes easy to access the abdominal cavity



# WOUND RETRACTOR

Laparoscopic Single Port



TYPE	PRODUCT NO.	WOUND DIAMETER	QUANTITY
String	WR-11-6060	60mm x 60mm	10ea/Box
	WR-11-9090	90mm x 90mm	10ea/Box
	WR-11-6090	60mm x 90mm	10ea/Box
Non-String	WR-12-6060	60mm x 60mm	10ea/Box
	WR-12-9090	90mm x 90mm	10ea/Box
	WR-12-6090	60mm x 90mm	10ea/Box
	WR-12-120	120mm x 120mm	5ea/Box
	WR-12-195	195mm x 195mm	5ea/Box