

A collection of white and blue laparoscopic surgical instruments, including trocars and graspers, arranged on a white background. The instruments are shown in various orientations, some in sharp focus and others blurred in the background. The text 'SEJONG MEDICAL' is overlaid in the upper right, and the company logo is in the bottom right.

**SEJONG MEDICAL**

*Laparoscopic Surgery System & Medical Equipment*


SEJONG MEDICAL CO.,LTD.  
**SEJONG  
MEDICAL**

# LAPORT-MU

Laparoscopic Trocar



TYPE	PRODUCT NO.	DIMENSION	QUANTITY
Optical	N21-0575	5mm x 75mm	10ea/Box
	N21-0510	5mm x 100mm	10ea/Box
	N21-1075	10mm x 75mm	10ea/Box
	N21-1010	10mm x 100mm	10ea/Box
	N21-1175	11mm x 75mm	10ea/Box
	N21-1110	11mm x 100mm	10ea/Box
	N21-1275	12mm x 75mm	10ea/Box
Blade	N21-1210	12mm x 100mm	10ea/Box
	N22-0575	5mm x 75mm	10ea/Box
	N22-0510	5mm x 100mm	10ea/Box
	N22-1010	10mm x 100mm	10ea/Box
	N22-1110	11mm x 100mm	10ea/Box
Bladeless	N22-1210	12mm x 100mm	10ea/Box
	N23-0575	5mm x 75mm	10ea/Box
	N23-0510	5mm x 100mm	10ea/Box
	N23-1010	10mm x 100mm	10ea/Box
	N23-1110	11mm x 100mm	10ea/Box
	N23-1210	12mm x 100mm	10ea/Box



# MINIPOINT

Laparoscopic Trocar



**Low Profile Design**  
Ideal for pediatric and petite patient

**Injury-Free Serration**  
Easy insert and prevents breakaway structure

**Needle Type**



Bladeless

TYPE	PRODUCT NO.	DIMENSION	QUANTITY
Non-Valve	TP52-0575	5mm x 75mm	10ea/Box
	TP52-0510	5mm x 100mm	10ea/Box
Valve	TP72-0575	5mm x 75mm	10ea/Box
	TP72-0510	5mm x 100mm	10ea/Box

# U-FLEX

Laparoscopic Trocar



**Smaller Fascial Defect**  
Soft cannula minimizes fascial defect

**Flexibility**  
Possible to insert curved instruments

**Needle Type**



Bladeless

**Low Profile Design**  
Ideal for pediatric and petite patient

TYPE	PRODUCT NO.	DIMENSION	QUANTITY
Non-Valve	T52-0575	5mm x 75mm	10ea/Box
Valve	T72-0575	5mm x 75mm	10ea/Box