




3rd Generation
**POWERED STAPLER
AND RELOADS**

**Freedom.
Control.
Confidence.**

Comprehensive reload range
optimises the balance between haemostasis and perfusion

Controlled and delicate powered jaw opening/closing
delivers maximum control and active tissue measurement

Built-in self-diagnostic
ensures readiness before every use



Teardrop groove guide
mitigates excessive tissue movement,
reducing the risk of malformed staples

Audible cues guide every step
offer clarity and confidence
to the entire surgical team

Tissue Sensing Technology
dynamically adjusts firing speeds to enhance
staple formation across varying tissue thicknesses

Improved Weight Distribution
provides greater control and precision,
reducing tension on
targeted and adjacent tissue

Thin/Medium	Thick/Extra Thick	Excessive
10 sec	20 sec	Stops

NEED MORE INFORMATION

Please contact us at
enquiries@lawmed.co.uk or **+44 (0)1932 260 838**.

Lawmed Ltd
Building 5 The Heights
Wellington Way
KT13 0NY
UK

www.lawmed.co.uk



3rd Generation
**POWERED STAPLER
AND RELOADS**

The Touchstone TSN Stapler brings innovation to your fingertips and delivers performance you can rely on.

Handle Code	Shaft Length (mm)	Shaft Diameter (mm)
TSN-S	65	12
TSN-M	160	12
TSN-L	260	12

Cartridge Code	Staple Line Length (mm)	Staple Height (mm)	Compatible Reload	Color
SRG-30W/45W/60W	30/45/60	2.0/2.5/3.0	SMG-30/45/60	○ White
SRG-30B/45B/60B	30/45/60	2.5/3.0/3.5	SMG-30/45/60	● Blue
SRG-45D/60D	45/60	3.0/3.5/3.8	SMG-45/60	● Gold
SRG-45G/60G	45/60	3.5/3.8/4.2	SMG-45/60	● Green

NEED MORE INFORMATION

Please contact us at
enquiries@lawmed.co.uk or **+44 (0)1932 260 838**.

Lawmed Ltd
Building 5 The Heights
Wellington Way
KT13 0NY
UK

www.lawmed.co.uk



3rd Generation
**POWERED STAPLER
AND RELOADS**

Reload Code	Staple Line Length (mm)	Staple Height (mm)	Color	Trocar Size
MG-30W/45W/60W	30/45/60	2.5	○ White	12
MG-30B/45B/60B	30/45/60	3.5	● Blue	12
MG-45D/60D	45/60	3.9	● Gold	12
MG-45G/60G	45/60	4.8	● Green	15
MT-45T/60T	45/60	5.0	● Black	15
FMG-30AR/ARH	30	2.0/2.5/3.0	● Tan	12
FMG-45AR/ARH	45	2.0/2.5/3.0	● Tan	12
FMG-60AR/ARH	60	2.0/2.5/3.0	● Tan	12
FMG-30PR/PRH	30	3.0/3.5/4.0	● Purple	12
FMG-45PR/PRH	45	3.0/3.5/4.0	● Purple	12
FMG-60PR/PRH	60	3.0/3.5/4.0	● Purple	12
FMG-45TR/60TR	45/60	4.0/4.5/5.0	● Black	15
SMG-30W/45W/60W	30/45/60	2.0/2.5/3.0	○ White	12
SMG-30B/45B/60B	30/45/60	2.0/2.5/3.0	● Blue	12
SMG-45D/60D	45/60	3.0/3.5/3.8	● Gold	12
SMG-45G/60G	45/60	3.5/3.8/4.2	● Green	12
SMT-45T/60T	45/60	4.2/4.5/5.0	● Black	15

NEED MORE INFORMATION

Please contact us at
enquiries@lawmed.co.uk or **+44 (0)1932 260 838**.

Lawmed Ltd
Building 5 The Heights
Wellington Way
KT13 0NY
UK

www.lawmed.co.uk