

Female Pelvic Organ Prolapse
Vaginal/Cervical Stump Prolapse




DynaMesh®-PRS soft and **DynaMesh®-PRS visible** implants are intended to be used as bridging material and reinforce the soft tissue of the vaginal walls as part of surgical treatment for apical pelvic organ prolapse.


DynaMesh®-PRS



Colpo-/cervicosacropexy
- unilateral



Colpo-/cervicosacropexy
- unilateral
- anterior/posterior mesh plasty (for concomitant cystocele/rectocele)



Colpo-/cervicosacropexy
- unilateral

DynaMesh®-PRS - Animation: Colposacropexy *https://de.dyna-mesh.com/Vi046xx*	
DynaMesh®-PRS - Animation: Colposacropexy *https://de.dyna-mesh.com/Vi048xx*	
DynaMesh® MRI - Animation: MRI Reconstruction with DynaMesh®-PRS visible *https://de.dyna-mesh.com/Vi067xx*	

hergestellt durch / manufactured by: FEG Textiltechnik Forschungs- und Entwicklungsgesellschaft mbH
Prager Ring 70, 52070 Aachen, Germany · www.dyna-mesh.com
made in Germany

CA_P1_PRS1_A_001_en001a_240319_d


Use and Properties

Product	DynaMesh®-PRS soft (1)	DynaMesh®-PRS visible (2)
Surgical Treatment	Apical Pelvic Organ Prolapse (Vaginal Stump / Cervical Stump)	
Surgical Approach	Minimally Invasive / Open	
Surgical Technique	Sacropexy	
Fixation	- Anterior longitudinal ligament: non-absorbable suture or tacks - Vaginal/cervical stump and vaginal walls: interrupted non-absorbable suture	
Smooth Warp-Knitted Selvedges	●	
Defined Elasticity	● [TR130]	
Visible Technology	○	●
Materials	- Polyvinylidene fluoride (PVDF) (CAS 24937-79-9) > 99% (w/w) (1) (2) - Triiron tetraoxide (CAS 1317-61-9) < 1% (w/w) (2)	
Polymer (Monofilament)	PVDF	
Biocompatibility	● [TR1]	
Ageing Resistance	● [2^{A}, 5^{VIT}, 27^{A}, 52^{VIT}, 93^{A}, 101]	
Effective Porosity	● High effective porosity reduces inflammation and the risk of excessive scar formation. [103^{P}, TR132]	
Klinge's Mesh Classification	Class 1a [102^{P}, TR132]	

Product Range

DynaMesh®-PRS soft	I) 02 cm x 16 cm	PV350216*F1*
	II) 03 cm x 23 cm	PV350323*F1*
DynaMesh®-PRS visible	III) 03 cm x 23 cm	PV750323*F1/F10*
	IV) 3.3 cm x 24 cm	PV750424*F1/F10*
	V) 04 cm x 20 cm	PV750420*F1/F10*
DynaMesh®-PRS soft	VI) 05 cm x 27 cm	PV350527*F1*

FX = X unit(s)/box (e.g. *F3* = 3 unit(s)/box)



● Applies to all product sizes
○ Does not apply
[#] Reference "#" (see "References")
[TR#] Internal test report (see "internal test report references")
Limitations "A" animal trial, "B" bench test, "VIT" in-vitro trial, "P" published results based on the analysis of human mesh explants, "PB" published results mainly based on bench tests

More information: *https://en.dyna-mesh.com/dynamesh-prs-gb*

Distributed by:
LAWMED
ADVANCING SURGICAL POSSIBILITIES TOGETHER
LAWMED Ltd.
Ground Floor, Building 5, The Heights
Wellington Way, Weybridge, KT13 0NY
United Kingdom
E: enquiries@lawmed.co.uk
T: +44 (0)1932 260 838
www.lawmed.co.uk

hergestellt durch / manufactured by: FEG Textiltechnik Forschungs- und Entwicklungsgesellschaft mbH
Prager Ring 70, 52070 Aachen, Germany · www.dyna-mesh.com
made in Germany

CA_P1_PRS1_A_002_en001a_250530_d